

# Fast and Efficient Deuteration in Continuous Flow Using the H-Cube<sup>®</sup> Advance Flow Reactor

## INTRODUCTION

Deuterium-labeled compounds are widely used as research tools in chemistry. Their importance lies in a number of applications, such as proving **reaction mechanisms**, the investigation of a compound's **pharmacokinetic properties**, internal standards in **mass spectrometry**, and compound structure determination in **NMR spectroscopy**.

Conventional techniques for the synthesis of deuterated compounds utilize D<sub>2</sub> gas as a deuterium source. However, there are drawbacks to utilizing deuterium gas on a laboratory scale, such as the **handling of the gas** itself. Other methods have been employed to overcome this difficulty, such as **catalytic H-D exchange** reactions between H<sub>2</sub> and D<sub>2</sub>O. However, these methods are time-consuming and do not produce high purity D<sub>2</sub>. They also require high pressure, the use of a special catalyst, or an excess amount of a strong base or acid<sup>1</sup>.

The **H-Cube<sup>®</sup> Advance** continuous flow system is capable of generating deuterium gas from the **electrolysis of D<sub>2</sub>O**, which is readily available in **99.98% purity** and is easy to handle.



## TECHNIQUE

High-purity deuterated compounds can be generated in high yield as long as a few guidelines are followed. Removing H<sub>2</sub>O from the water reservoir and the system is crucial. We recommend flushing in Manual operation mode.

**Step 1:** Clean and dry the inside of the water reservoir with a paper towel.

**Step 2:** Fill the water reservoir with D<sub>2</sub>O. Place the free end of the water inlet tube into the reservoir.

**Step 3:** Select internal H<sub>2</sub> on the device's screen, and wait for the pressure in the electrolytic cell to build up.

**Step 4:** Prime the HPLC pump and both the solvent and reagent inlet tubes using a syringe adapter.

**Step 5:** Place the outlets into a product and waste vessel respectively.

**Step 6:** Tap the waste collection button (🗑️) on the screen. Start the HPLC pump at a set flow rate to wash the system with clean solvent.

**Step 7:** Set the desired gas flow rate and start the gas flow. Let the system run for approximately 10-15 minutes.

**Step 8:** Empty the water reservoir and refill it with fresh D<sub>2</sub>O, making sure the system is completely purged of trace H<sub>2</sub>O from previous runs.

After the system has been properly flushed, the reaction can be started.

## RECOMMENDATIONS

- Only use an **aprotic solvent** to avoid H-D exchange.
- Do not use hydrogen saturated catalysts such as **Raney catalysts** for the reaction.


## LITERATURE

There are two literature examples of where the H-Cube® Advance has been utilized for the deuteration of compounds.

Fülöp et al. performed deuteration on a series of unsaturated compounds<sup>1</sup>.

Optimum conditions were first explored using **cinnamic acid** as a standard.

Early reactions using methanol and **10% Pd/C** led to only a **30% incorporation** of deuterium into the molecule. Changing from methanol to an aprotic solvent, **ethyl acetate**, increased **deuterium incorporation to 70%**.



The catalyst was then changed to a less active catalyst (**5% Pd/BaSO<sub>4</sub>**) to reduce the deuteration of the phenyl ring and selectively deuterate the double bond. This catalyst change increased **deuterium content to 95%**. The reaction was carried out at room temperature. Pressures in the 40-60 atm. range and flow rates between 0.7-2 mL/min did not affect yield or deuterium levels.

Once the optimized conditions had been found, the same conditions were applied to a series of other unsaturated compounds. The results are displayed in Table 1.

As you can see from the results, the products were synthesized in near quantitative yield with a high deuterium incorporation. No purification was necessary. The structures included **foldamer building blocks**, so there is potential for structure elucidation via deuteration where bacteria labeling is not possible. The D<sub>2</sub>O consumption was very low (4.41 µL/min), which is a much higher deuterium efficiency when compared to other methods.

The other paper is from Kappe et al. and describes the deuteration of **ethyl cinnamate**<sup>2</sup>. Using a flow rate of 1 mL/min, 10% Pd/C, room temperature, and Full H<sub>2</sub> mode, the product was successfully synthesized in **92% yield** with a **95%** degree of **deuterization**.

Entry	Substrate	Product	D <sub>2</sub> <sup>a</sup> (%)	Yield <sup>b</sup> (%)	Entry	Substrate	Product	D <sub>2</sub> <sup>a</sup> (%)	Yield <sup>b</sup> (%)
1			99	99	5			96	99
2			97	98	6			95	99
3			93	97	7			97	95
4			96	98	8			98	98

<sup>a</sup> Deuterium content in %. <sup>b</sup> Isolated yield.

## CONCLUSION

Utilizing D<sub>2</sub>O instead of H<sub>2</sub>O, the H-Cube® Advance is able to deuterate compounds in **high yield** and with a **deuteration incorporation of >95%**. These results offer up the H-Cube® Advance as a reliable alternative to other deuteration methods.

## REFERENCES

[1] Mándity, I.M.; Martinek T.A.; Darvas, F.; Fülöp, F.; *Tet. Lett.*; **2009**; 50; 4372-4374.

[2] Irfan, M.; Petricci, E.; Glasnov, T.; Taddei, M.; Kappe, O.; *Eur. J. Org. Chem.*; **2009**; 1327-1334.

For more information, please visit  
[www.thalesnano.com](http://www.thalesnano.com)



ThalesNano Inc.  
Záhony utca 7.  
H-1031 Budapest | Hungary  
Email: [sales@thalesnano.com](mailto:sales@thalesnano.com)  
[www.thalesnano.com](http://www.thalesnano.com)

US Office  
50 S. Penn St. Suite B-2  
Hatboro PA. 19040 USA  
Phone: 215-534-3365  
E-mail: [USAsales@thalesnano.com](mailto:USAsales@thalesnano.com)

