



BACK PRESSURE MODULE 300

Back Pressure Module 300™

Widening the available
chemical space



- Pressure up to 300 bar (4650 psi)
- Improved stability
- Outlet switch valve



ABOUT THE SYSTEM

With our new back pressure regulator, reactions can be performed up to 300 bar (4650 psi) pressure in a safe and controlled manner. Improved stability is ensured by the two separate system pressure valves.

The Back Pressure Module 300™ is an easy-to-use motorized back pressure regulator designed to control the pressure of the connected flow chemical system in a reliable and convenient manner. It is equipped with a touch screen on which customers can set up the desired conditions.

FEATURES

Working pressure range	60–300 bar (870–4650 psi)
Applicable flow rates:	0.3–10 mL/min for liquids
Pressure precision:	± 2 bar (30 psi)



For more information, please visit
www.thalesnano.com
Linkedin: /company/thalesnano-inc-
Twitter: /thales_nano
Instagram: /thalesnano_inc/
Facebook: /ThalesNano/

ThalesNano Inc.
Zahony u. 7.
H-1031 Budapest
Hungary
Tel.: +36 1 880 8500
Fax: +36 1 880 8501
E-mail: sales@thalesnano.com

US Office
50 S. Penn St. Suite B-2
Hatboro
PA. 19040
USA
Phone: 215-534-3365
E-mail: USAsales@thalesnano.com